



FEGURA SIL CLEAR

Transparent duplicating silicone, shore A 25

For the processing of light-curing dental acrylics, addition-curing

Fegura Sil clear can be used to produce silicone moulds for the processing of light-curing dental acrylics (e.g. composite materials). Light hardening is performed through the silicone mould which prevents the formation of a dispersion layer. Thus, the surface becomes tighter and can be polished more easily.

- processing with duplicating flask
- also for demonstration models

1 kg A + 1 kg B (REF 1095)



FEGURA SIL GLASS

Transparent cartridge silicone

For the processing of light-curing dental acrylics, addition-curing

Fegura Sil GLASS can be used to produce silicone moulds for the processing of light-curing dental acrylics (e.g. composite materials). Light hardening is performed through the silicone mould which prevents the formation of a dispersion layer. Thus, the surface becomes tighter and can be polished more easily.

Fegura Sil GLASS is available in two hardness grades:

- Shore A 55 - regular
- Shore A 75 - super hard

Set: 2 x 50 ml, 12 mixing cannulas

2 x 50 ml set 55 regular (REF FSG0120)

2 x 50 ml set 75 super hard (REF FSG0100)

