



FEGURA SIL GINGI EXTRA HARD

Gingival mask is hard and allows for easy milling

For producing gingival masks for dental implants

Is directly injected into the impression. Easy milling of the material is possible thanks to the high hardness, thus achieving precise definition of margins.

- no costly and time-consuming corrections on completed restorations
- reproduction of the natural gingival situation
- easy and economical to use

Set: 2 x 50 ml, mixing cannulas and intra-oral tips (12 pieces each), separating agent

Refill pack: 2 x 50 ml, mixing cannulas and intra-oral tips (12 pieces each)

2 x 50 ml set (REF 2015)

2 x 50 ml refill pack (REF 2016)

15 ml separating agent (REF 2012)

FEGURA SIL GINGI

Highly elastic gingival mask

For producing gingival masks for saw-cut models

Has the same colour, hardness and flexibility as natural dental tissue. Prosthetic restorations can be perfectly adapted due to the detailed reproduction of the natural gingival situation, especially concerning the functional design of interdental areas and the correct placement of attachments.

The silicone impression, into which Fegura Sil GINGI is injected, can be made from the kneading silicone Fegura Sil putty.

- no costly and time-consuming corrections on completed restorations
- reproduction of the natural gingival situation
- easy and economical to use

Set: 2 x 50 ml, mixing cannulas and intra-oral tips (12 pieces each), separating agent

Refill pack: 2 x 50 ml, mixing cannulas and intra-oral tips (12 pieces each)

2 x 50 ml set (REF 2010)

2 x 50 ml refill pack (REF 2011)

15 ml separating agent (REF 2012)