



Website



YouTube



Facebook




Changzhou Sifary Medcial Technology Co., Ltd
Add: NO.99, Qingyang Road, Xuejia County, Xinbei
District, Changzhou City, 213000 Jiangsu, China
Email: ivy@sifary.com
Web: www.eighteeth.com



Landlink GmbH

Add: Dorfstrasse 2/4, 79312, Emmendingen, Germany
Email: info@landlink.eu

DISTRIBUTED BY



TIP BOOK



CATALOG

Scaling

G1

Order Code: 6680523



G1-S

Order Code: 6680617



G2

Order Code: 6680524



G2-S

Order Code: 6680618



G4

Order Code: 6680526



G4-S

Order Code: 6680619



G5

Order Code: 6680527



G5-S

Order Code: 6680620



Periodontic

P1

Order Code: 6680528



5

P1-S

Order Code: 6680621



P3

Order Code: 6680529



6

P3-S

Order Code: 6680622



P4

Order Code: 6680530



7

P4-S

Order Code: 6680623



7

P52

Order Code: 6680533



8

P52-S

Order Code: 6680626



P50L

Order Code: 6680531



9

P50L-S

Order Code: 6680624



Periodontic

P50R

Order Code: 6680532



10

P50R-S

Order Code: 6680625



P53L

Order Code: 6680534



11

P53L-S

Order Code: 6680627



P53R

Order Code: 6680535



12

P53R-S

Order Code: 6680628



12

P59

Order Code: 6680536



13

P59-S

Order Code: 6680630



P56

Order Code: 6680544



14

P56-S

Order Code: 6680629



Periodontic

P3D

Order Code: 6680654



15

P3-SD

Order Code: 6680636



P4D

Order Code: 6680631



16

P4-SD

Order Code: 6680637



Endodontic

E1

Order Code: 6680537



17

E1-S

Order Code: 6680656



E3

Order Code: 6680538



18

E3-S

Order Code: 6680657



E4

Order Code: 6680539



19

Endodontic

E4-S

Order Code: 6680658



19

E5

Order Code: 6680540



20

E5-S

Order Code: 6680659



E7

Order Code: 6680542



21

E7-S

Order Code: 6680660



E8

Order Code: 6680608



22

E8-S

Order Code: 6680614



E9

Order Code: 6680609



23

E9-S

Order Code: 6680615



E10

Order Code: 6680610



24

Periodontic

E10-S

Order Code: 6680616



24

E3D

Order Code: 6680651



25

E3-SD

Order Code: 6680638



E30

Order Code: 6680611



26

E30-S

Order Code: 6680661



Endodontic

E31

Order Code: 6680612



27

E31-S

Order Code: 6680662



E4D

Order Code: 6680652



28

E4-SD

Order Code: 6680639



E5D

Order Code: 6680653



29

Endodontic

E5-SD

Order Code: 6680640



29

E6D

Order Code: 6680632



30

E6-SD

Order Code: 6680641



E7D

Order Code: 6680633



31

E7-SD

Order Code: 6680642



E8D

Order Code: 6680634



32

E8-SD

Order Code: 6680643



Cavity Preparation

TC1D

Order Code: 6680648



33

TC1-SD

Order Code: 6680650










SCALING

EMS® **G1**
Satelec® **G1-S**

Order Code **6680523**
Order Code **6680617**


Function Description: All-round removal of supragingival calculus



Mode: G
Power: 1-10



Recommended
Service Life
200 Teeth








Water Flow Rate
mid&high

EMS® **G2**
Satelec® **G2-S**

Order Code **6680524**
Order Code **6680618**


Function Description: Removal of significant supragingival calculus



Mode: G
Power: 1-10



Recommended
Service Life
200 Teeth








Water Flow Rate
mid&high

EMS® **G4**
 Satelec® **G4-S**

Order Code **6680526**
 Order Code **6680619**






Function Description: Used for supragingival and interdental space scaling



EMS® **G5**
 Satelec® **G5-S**

Order Code **6680527**
 Order Code **6680620**

Function Description: Removal of significant supragingival calculus










PERIODONTIC

EMS® **P1**
Satelec® **P1-S**

Order Code **6680528**
Order Code **6680621**

Function Description: Removal of subgingival calculus and periodontal irrigation



P
Mode: P
Power: 1-10

Recommended
Service Life
150 Teeth

Water Flow Rate
mid&high






EMS® P3

Satelec® P3-S

Order Code 6680529

Order Code 6680622


Function Description: Removal of calculus deep in periodontal pocket and its irrigation



Mode: P
Power: 1-6



Recommended
Service Life
150 Teeth



Water Flow Rate
mid&high






EMS® P4

Satelec® P4-S

Order Code 6680530

Order Code 6680623


Function Description: Removal of calculus in shallow pockets



Mode: P
Power: 1-6



Recommended
Service Life
200 Teeth





Water Flow Rate
mid&high

EMS® P52
Satelec® P52-S

Order Code 6680533
Order Code 6680626



Function Description: Supragingival scaling



EMS® P50L
Satelec® P50L-S

Order Code 6680531
Order Code 6680624


Function Description: Used for subgingival calculus probing, subgingival proximal surface scaling and root planing



EMS® P50R
Satelec® P50R-S

Order Code **6680532**
 Order Code **6680625**

Function Description: Used for subgingival calculus probing , subgingival proximal surface scaling and root planing



P
 Mode: P
 Power: 1-10


Recommended Service Life
 150 Teeth

Water Flow Rate
 high

EMS® P53L
Satelec® P53L-S

Order Code **6680534**
 Order Code **6680627**

Function Description: Used for subgingival calculus probing , subgingival proximal surface scaling and root planing



P
 Mode: P
 Power: 1-10


Recommended Service Life
 150 Teeth

Water Flow Rate
 mid&high

EMS® P53R
Satelec® P53R-S

Order Code **6680535**
 Order Code **6680628**


Function Description: Used for subgingival calculus probing , subgingival proximal surface scaling and root planing



Mode: P
 Power: 1-10



Recommended
 Service Life
 150 Teeth




Water Flow Rate
 mid&high

EMS® P59
Satelec® P59-S

Order Code **6680536**
 Order Code **6680630**


Function Description: Used for subgingival calculus probing and scaling



Mode: P
 Power: 1-6



Recommended
 Service Life
 150 Teeth




Water Flow Rate
 high

EMS® **P56**
 Satelec® **P56-S**

Order Code **6680544**
 Order Code **6680629**

Function Description: Removal of significant supragingival calculus and plaque



P
 Mode: P
 Power: 1-10


Recommended Service Life
 80 Teeth

Water Flow Rate
 high

EMS® **P3D**
 Satelec® **P3-SD**

Order Code **6680654**
 Order Code **6680636**

Function Description: Used to level off the surface of endo after the periodontal flap surgery



P
 Mode: P
 Power: 1-6

Recommended Service Life
 50 Teeth

Water Flow Rate
 mid&high


EMS® P4D

Order Code 6680631

Satelec® P4-SD

Order Code 6680637

Function Description: For locating root canal and removal of calcification on the one-third of crown



Mode: P
Power: 1-6



Recommended
Service Life
50 Teeth



Water Flow Rate
low, mid & high

ENDODONTIC


EMS® E1

Satelec® E1-S

Order Code 6680537

Order Code 6680656


Function Description: 120°angle holder for 0.7mm-0.8mm endo file for root canal cleaning



Mode: E
Power: 1-3



Recommended
Service Life
500 Teeth



Water Flow Rate
mid&high


EMS® E3

Satelec® E3-S

Order Code 6680538

Order Code 6680657


Function Description: For lateral condensation of gutta percha by heating effect



Mode: E
Power: 1-3



Recommended
Service Life
500 root canals



Water Flow Rate
No water


EMS® E4

Satelec® E4-S

Order Code 6680539

Order Code 6680658


Function Description: Removal of foreign filling material during re-treatment



Mode: E
Power: 1-6



Recommended
Service Life
100 root canals



Water Flow Rate
No water


EMS® E5

Satelec® E5-S

Order Code 6680540

Order Code 6680659


Function Description: Removal of foreign filling material during re-treatment



Mode: E
Power: 1-6



Recommended
Service Life
100 root canals



Water Flow Rate
No water


EMS® E7

Satelec® E7-S

Order Code 6680542

Order Code 6680660


Function Description: Removal of broken instrument



Mode: E
Power: 1-3



Recommended
Service Life
50 root canals



Water Flow Rate
mid&high


EMS® E8

Satelec® E8-S

Order Code 6680608

Order Code 6680614


Function Description: Used to brush the canal walls during re-treatment



Mode: E
Power: 1-6



Recommended
Service Life
50 root canals



Water Flow Rate
low,mid&high


EMS® E9

Satelec® E9-S

Order Code 6680609

Order Code 6680615


Function Description: Removal of broken instrument



Mode: E
Power: 1-6



Recommended
Service Life
500 root canals



Water Flow Rate
No water


EMS® E10

Satelec® E10-S

Order Code 6680610

Order Code 6680616

Function Description: Used for irrigation after preparation



Mode: E
Power: 1-3



Recommended
Service Life
50 root canals



Water Flow Rate
mid&high


EMS® E30

Satelec® E30-S

Order Code 6680611

Order Code 6680661


Function Description: Removal of deep foreign particles in canal



Mode: E
Power: 1-3



Recommended
Service Life
50 root canals



Water Flow Rate
low&mid


EMS® E31

Satelec® E31-S

Order Code 6680612

Order Code 6680662


Function Description: Removal of foreign particles in one-third of the canal



Mode: E
Power: 1-3



Recommended
Service Life
50 root canals




Water Flow Rate
low&mid

EMS® E3D
Satelec® E3-SD

Order Code 6680651
Order Code 6680638

Function Description: Enlargement of the canal walls to remove the broken instrument



Mode: E
Power: 1-3



Recommended
Service Life
50 root canals




Water Flow Rate
low, mid & high

EMS® E4D
Satelec® E4-SD

Order Code 6680652
Order Code 6680639


Function Description: Removal of hard material in canal during re-treatment



Mode: E
Power: 1-3



Recommended
Service Life
50 root canals




Water Flow Rate
No water

EMS[®] E5D
Satelec[®] E5-SD

Order Code 6680653
Order Code 6680640


Function Description: Removal of hard material in canal during re-treatment



Mode: E
Power: 1-3



Recommended
Service Life
50 root canals




Water Flow Rate
No water

EMS[®] E6D
Satelec[®] E6-SD

Order Code 6680632
Order Code 6680641

Function Description: Removal of foreign material in cavity to access canal orifice



Mode: E
Power: 1-10



Recommended
Service Life
50 root canals




Water Flow Rate
mid&high

EMS® E7D
Satelec® E7-SD

Order Code 6680633
Order Code 6680642

Function Description: Removal of calcification and filling material in cavity to access canal orifice during re-treatment



Mode: E
Power: 1-3



Recommended
Service Life
30 root canals




Water Flow Rate
mid&high

EMS® E8D
Satelec® E8-SD

Order Code 6680634
Order Code 6680643


Function Description: Used for root canal retrogression, efficient root apical polishing



Mode: E
Power: 1-6



Recommended
Service Life
50 root canals



Water Flow Rate
low,mid&high

CAVITY PREPARATION


EMS® TC1D

Satelec® TC1-SD





Order Code 6680648

Order Code 6680650

Function Description: Removal cervical caries and caries on the occlusal surface



Mode: P
Power: 1-10



Recommended
Service Life
30 Teeth



Water Flow Rate
low, mid & high