



FEGU-SOLV

powerful plaster removing agent

To remove any plaster residues

- cleaning impression trays, mixing beakers and acrylic dentures
- control of the activity by means of bubbling
- especially effective in the ultrasonic bath

Technical data	
Colour	green

1000 ml (REF 5095)
5000 ml (REF 5097)

FEGU-SOLV SUPER

Very strong plaster solvent,

To remove any plaster residues

- increased dissolving capacity
- for impression trays, mixing beakers and acrylic dentures
- control of the activity by means of bubbling
- especially effective in the ultrasonic bath

Technical data	
Colour	yellow

1000 ml (REF 5092)
5000 ml (REF 5093)